

Water Bath
BM42/82



BM42/82

- Simple water bath with stick thermometer
- Beaker can be used up to 2L (BM42), 5L (BM42)

■ Specifications

Model	BM42	BM82
Standard		
Operating temperature range	Room temp. +5deg.C to Boiling temp.	
Temp adjustment accuracy	±2.5 deg.C	
Temp. setting	Analog set-up	
Heater	500W	1kW
Temp. control system	Bimetal temp. control + Electronical control	
Bath capacity/Material	3.7L/SUS304	12.3L/SUS304
Internal dimension	Dia 200 x H120 (mm)	Dia 280 x H200 (mm)
External dimension	W202 x D250 x H120(mm)	W315 x D360 x H200(mm)
Power	AC100V, 50/60Hz, 5A	AC100V, 50/60Hz, 10A
Weight	Approx. 1.5kg	Approx. 2.6kg
Accessory	Stick thermometer x 1pc.	

The performance shown here is the value under the condition of the ambient temperature of 23deg.C±5deg.C, humidity of 65%RH±20%.

Model	
42	82
Operating temp. range	
Room Temp. +5 to Boiling Temp.	
Temp adjustment accuracy	
±2.5 deg.C	
Bath Capacity	
3.7L	12.3L