

Shaking Unit
BN300



- This unit can be operated with Water Bath (Model: BQ, BV, BK, BA or the other baths if the capacity is suitable) for shaking purposes with easy put-on and take-off.
- Loaded with “Shaking Width Changeable Function” with digital setting and display.

Shaking Width / Times
10 to 30mm changeable / 20 to 160 times /min.

■ Specifications

Model	BN300
Performance	
Shaking Times	20 to 160 times/min.
Shaking Width	10 to 30mm changeable
Setting and Display	Digital
Shaking Motor	Condenser Motor, 25W
Safety Device	
Safety Countermeasures	Motor Overloading Protection, Power Fuse
Standard	
Number of Sample	Test Tube (15mm dia. x 80pcs.) Flask (300ml x 6pcs.) Flask (500ml x 3pcs.)
Spring Pitch	15mm pitch (square)
External Dimensions	W340 x D644 x H250mm
Power	AC100V, 50/60Hz, 1A (with no breaker)
Weight	Approx. 10kg



Combination Usage: BN300 + BV300



BN300 +BK400



BN300+BF401+BZ300