

# Shaking Unit

## BN300

### Shaking Width / Times

10 to 30mm changeable / 20 to 160 times /min.



- This unit can be operated with Water Bath (Model: BQ, BV, BK, BA or the other baths if the capacity is suitable) for shaking purposes with easy put-on and take-off.
- Loaded with “Shaking Width Changeable Function” with digital setting and display.

### ■ Specifications

Model	BN300
<b>Performance</b>	
Shaking Times	20 to 160 times/min.
Shaking Width	10 to 30mm changeable
Setting and Display	Digital
Shaking Motor	Condenser Motor, 25W
<b>Safety Device</b>	
Safety Countermeasures	Motor Overloading Protection, Power Fuse
<b>Standard</b>	
Number of Sample	Test Tube (15mm dia. x 80pcs.) Flask (300ml x 6pcs.) Flask (500ml x 3pcs.)
Spring Pitch	15mm pitch (square)
External Dimensions	W340 x D644 x H250mm
Power	AC100V, 50/60Hz, 1A (with no breaker)
Weight	Approx. 10kg



Combination Usage: BN300 + BV300



BN300 +BK400



BN300+BF400+BZ300