

BO300 has been discontinued.



- Digital Temp. Setting and Display.
- Loaded with the more durable heater made of SUS.

Model
BO300
Operating Temperature Range
50deg.C to +200deg.C

■ Specifications

Model	BO300
Performance	
Operating Temp. Range	50 to 200deg.C
Temp. Adjustment Accuracy	±0.5deg.C(at 200deg.C)
Temp. Distribution Accuracy	±7deg.C (at 200deg.C)
Time to Reach	20 to 200deg.C Approx. 70min.
Components	
Temp. Control System	Micro Computer PID Control
Heater	SUS316L, 1.2kW
Safety Device	
Safety Countermeasures	Self Diagnosis Function, Electric Leakage Breaker, Double Heating Preventor, Empty Boiling Preventor
Standard	
Effective Dimension	Φ45, 75, 105, 135, 195, 225mm
Inner Dimensions	W270 x D270 x H170mm
External Dimensions	W340 x D440 x H310mm
Inner Capacity	Approx. 12L
Power	AC100V, 50/60Hz, 12A
Weight	Approx. 10kg
Accessory	Shelf 1pc., Ring-Shape Lid 1set, Pole 1pc., Clamp 1pc., Ground Adapter 1pc.