

Oil Bath**BO300**

Model
BO300
Operating Temperature Range
50deg.C to +200deg.C

FEATURES

- Digital Temp. Setting and Display.
- Loaded with the more durable heater made of SUS.

**Specifications**

Model	BO300
■ Performance	
Operating Temp. Range	50 to 200deg.C
Temp. Adjustment Accuracy	±0.5deg.C(at 200deg.C)
Temp. Distribution Accuracy	±7deg.C (at 200deg.C)
Time to Reach	20 to 200deg.C Approx. 70min.
■ Components	
Temp. Control System	Micro Computer PID Control
Heater	SUS316L, 1.2kW
■ Safety Device	
Safety Countermeasures	Self Diagnosis Function, Electric Leakage Breaker, Double Heating Preventor, Empty Boiling Preventor
■ Standard	
Effective Dimension	Φ45, 75, 105, 135, 195, 225mm
Inner Dimensions	W270 x D270 x H170mm
External Dimensions	W340 x D440 x H310mm
Inner Capacity	Approx. 12L
Power	AC100V, 50/60Hz, 12A
Weight	Approx. 10kg
Accessory	Shelf 1pc., Ring-Shape Lid 1set, Pole 1pc., Clamp 1pc., Ground Adapter 1pc.