


Constant Temperature Water Bath

BS200/400/600/660

Model			
200	400	600	660
Operating Temperature Range			
Room Temp. +5 to Boiling Temp. (Water)			
Temp. Distribution Accuracy			
±3deg.C (at 70deg.C)			

FEATURES

- BS200: with 4 kinds of Baskets. Max. 5L Beaker can be set by removing Lid.
- BS600: Can use 6pcs. of Condensers simultaneously.
- BS400: Can use 4pcs. of Condensers simultaneously.
- BS660: Can use 6pcs. of Big/Small Container and can adjust the sink level.



BS200/400/660/600

Control Panel




Useable for various containers



Removable Lid



Easy to change Water Level by adjusting Over Flow



Well-Holding Clamp

Specifications

Model	BS200	BS400	BS600	BS660
■ Performance				
Operating Temp. Range	Room Temp. +5deg.C to Boiling Temp. (Water)			
Temp. Adjustment Accuracy	±3deg.C (at 70deg.C)			
Temp. Distribution Accuracy	±3deg.C (at 70deg.C)			
Time to Reach to Constant Max. Temp.	Approx. 35min.		Approx. 60min.	Approx. 50min.
■ Component				
Inner Bath	Stainless Steel			
Heater	Copper Pipe Heater (Nickel Plated)			
	900W	1.4kW		1.3kW
■ Safety				
Safety Countermeasure	Overheat Prevention System (Liquid Expansion Type, Self Sustaining Relay, Alarm Lamp On) Electric Leakage Breaker with Over Current Protection (Power Lamp On at the time System On)			
■ Standard				
Internal Dimensions	φ200 x D150mm	W328 x D302 x H101mm	W790 x D150 x H102mm	W502 x D302 x H101mm
Effective I.D.	φ186 x D119mm	W298 x D260 x H73mm	W776 x D140 x H74mm	W468 x D260 x H73mm
External Dimension	W300 x D364 x H217mm	W390 x D438 x H214mm	W864 x D288 x H215mm	W565 x D437 x H214mm
Bath Capacity	Approx. 4.7L	Approx. 9L	Approx. 11L	Approx. 14L
Power (50/60Hz)	AC100V/9A	AC100V/14A		AC100V/13A
Opening Size, Q'ty	4 pcs. with baskets	φ115mm, 4pcs	φ115mm, 6pcs	
Weight	Approx. 5kg	Approx. 9.5kg	Approx. 16kg	Approx. 12kg
Accessory	Ring (Not including BS200), Lid, Container Fixing Clamp, Stand, Bar for Clamp Fixing (Not including BS200), Connector for Stand, Thermometer (0 to +100deg.C Alcohol Type), Bottom Plate, Basket (Only for BS200)			

■ Optional items

Description	Model No.	Product Code
Cooling Water Distribution System for BS600, Supply/Drain Tubeφ16mm, Cock φ8.5mm	BS Series	222086
Cooling Water Distribution System for BS400, Supply/Drain Tubeφ12mm, Cock φ8.5mm	BS Series	222081