

## Fume Hoods

### Model KFS



Non-standard model with improved safety and working efficiency.

- The width and work surface can be selected according to the specific application.  
Width: 1,200, 1,500, 1,800mm  
Work surface: ceramitite, epoxy, and ceramic steel plate
- Optional air quantity controller, wind speed meter or air foil may be installed.

### Achieved safety, amenity and energy-savings Combining with an air quantity controller and a KF converter



The air quantity controller has two control modes: fixed wind speed and fixed air quantity. Output signals from one of these modes will be processed in the KF converter to be used for controlling the central exhaust or for monitoring the status of several KFS units from the central monitor room.

