

PCV-845/1305/1605/1915 have been discontinued.



Energy saving is achieved by adopting new type turbo fan. It minimized the intrusion of dust with vacuum structure around the HEPA filter.

Standard Type

■ Specifications

Product Code		100000	100000	100000	100000	100000	100000	100000	100000	
Model		PCV	845BN	845CN	1305BN	1305CN	1605BN	1605CN	1915BN	1915CN
Dimension (mm)	Height	H2	1,860	2,015	1,860	2,015	1,860	2,015	1,860	2,015
	Width	W1	840		1,300		1,600		1,910	
Performance	Wind Velocity	Initial: Approx. 0.4 m/s Telophase: Approx. 0.2 m/s								
	Cleanliness	Class 10 (Particle Dia.: 0.5µm at surrounding class 100,000)								
	Dust Collecting Efficiency	99.99% or more (with Particle Dia of 0.3µm)								
	Noise [dB (A) 50/60 Hz]	53/54 (with PCJ-1305BN)								
	Power Source	100VAC Single Phase 50/60 Hz								
	Consumption Power (50/60 Hz)	135/150W		180/210W		245/280W		255/295W		
	Weight (kg)	150	160	190	200	220	230	300	310	