

Low Constant Temperature Water Bath

BU200/300

Model	
200	300
Operating Temp. Range	
-10 to +80deg.C	-20 to +80deg.C
Inner Capacity	
13L	30L
Type	
Desk Top	Floor Stand

FEATURES

- Constant Low Temp. Operation by Combining with “Thermomate”



Specifications

Model	BU200	BU300
■ Performance		
Operating Temp. Range	-10 to +80deg.C	-20 to +80deg.C
Cooling Capacity (J/S)	350 (300Kcal/h) at 20deg.C	465 (400Kcal/h) at 20deg.C
	290 (250Kcal/h) at 10deg.C	400 (350Kcal/h) at 0deg.C
Time to Reach	Approx. 90min. (20deg.C -> -10deg.C)	Approx. 150min. (20deg.C -> -20deg.C)
■ Component		
Refrigerator	Cooling Air Total Closed Type150W	Cooling Air Total Closed Type400W
Coolant	R22	
Others	Service ReceptacleAC100V12A(for Thermomate BF Model)	N/A
	Drain	
■ Safety Device		
Safety Countermeasure	Electric Leakage Breaker, Refrigerator Over Load Relay	
■ Standard		
Internal Dimension	W157 x D220 x H136mm	W229 x D360 x H184mm
External Dimension	W653 x D478 x H314mm	W385 x D575 x H839mm
Water Bath Capacity	13L	30L
Power (Set with BF Motel)	AC100V15A	AC100V6A
Weight	Approx. 39kg	Approx. 52kg
Accessory	Nil	Bath Cover for Low Temperature

■ Optional Item (for BU200)	Model
Bath Cover for Low Temp.	OBU10
Stand (with Height Adjustment)	OBW10