

Loop Cinerator

SL-21

Model

SL-21

Features

- Pt. Ear, Needle is sterilized in Combustion Tube
- Safety Operation by Heater Sterilization
- Only 1/6 Heat Amount comparing to Gas Burner
- Good Heat-resistance
- Loaded with Quartz Glass which evenly Heatable
- Efficient for big volume and for Continuous Sterilizing



Flexible Angle Changing



Pt Ear can be put in holder holes

Consumables



Specifications

Model	SL-21
■ Standard	
Internal Sterilizing Temp.	800 to 850 deg.C
Time to Reach	10min. (to Sterilizing Temp.)
Cinerating Port Diameter	15mm
Heater	130W
Power	AC100V, 50/60Hz, 1.5A
External Dimensions.	W95 x D170 x H200mm
Weight	1.5kg