


Mill**UB32**

Model
UB32
Roller Revolution
200 to 600 rpm

Features

- Automatic Operation by Timer
- Easy Changeable Pot (9 to 30cm) to meet Sample Size
- Roller Revolution is controllable by the knob
- Stable Rolling by Flexible Pitch Adjustment based on the Pot Size and Special Rubber Roller

**Specifications**

Model	UB32
■ Standard	
Roller	Natural Rubber Made, $\Phi 50$ x Length 650mm
Roller Revolution	Approx. 200 to 600rpm, Step-less Revolution Control
Motor	Condenser Type, 4 poles 200W, Feedback Control
Timer	60/72min. (50/60Hz)
External Dimensions	W1,010 x D410 x H380mm
Power	AC100V, 3.4/3.7A (50/60Hz)
Weight	Approx. 51kg

■ Optional Item

Porcelain Pot (with Porcelain Pole)	Qty.	Revolution	Dimensions
	4	111 to 333rpm	90mm
	4	83 to 250rpm	120mm
	3	67 to 200rpm	150mm
	2	58 to 167rpm	180mm
	2	47 to 143rpm	210mm
	1	42 to 125rpm	240mm
	1	33 to 100rpm	300mm