

Mill**UT22**

Model	UT22
Pestle Revolution	50 to 110 rpm
Mortar External Diameter	Φ200mm

Features

- Continuous and Changeable Revolution
- With Timer
- Speedy Milling by opposite rotation of Two Pestles and a mortar
- With Dust Prevention Cover

**Specifications**

Model	UT22
■ Standard	
Revolution (50/60Hz)	Pestle 55 to 110rpm Continuous Variable (Motor: Inverter Type) Mortar 6/7rpm
Mortar	External 200mm Porcelain Made
Pestle	Φ30mm, Alumina Made, 2pcs.
Timer	60min. (Spring Type)
External Dimensions	W372 x D444 x H800mm
Power	AC100V, 50/60Hz, 1A
Weight	Approx. 31kg
Accessory	Porcelain Mortar 1pc., Spatula 1pc. Alumina Pestle 2pcs.

■ Optional Item
Alumina Made Pestle
Alumina Made Mortar (O.D. 200mm)
Porcelain Made Mortar
Agate Made Mortar (I.D. 140mm)
Agate Made Pestle (1pc.)
Agate Made Mortar Fixing Holder
Spatula (Not including Agate Mortar)