

Water Purifier

WB-21

Model	WB-21
Distillation Capacity	1.8 L/h
Pure Water Collection Type	Distillation Type

Features

- Design is conforming to JIS K0101 and K0102
- Suitable for Water Quality Test as the Purified Water for Mixing
- Raw Water Tank/Pure Water Tank is detachable
- Initial Produced Water can be automatically drained to remove Volatile Contamination



Specifications

Model	WB-21
■ Performance	
Collection Type	Distillation
Collection Pure Water	Distilled Water
Distilled Water Production Amount	Approx. 1.8L/h
Distilled Water Supply	Approx. 3L/min. (with Full Tank)
■ Components	
Boiler / Condenser	Super Hard Glass
Heater	Built-in Quartz Glass Tube, Bobbin Coil Heater
Distilled Water Tank	Polyethylene Made 10L (with 0.2μm Bacteria Elimination Filter)
Raw Water Tank	Polyethylene Made 10L
Power	AC100V, 50/60Hz, 15A
External Dimensions	W600 x D550 x H1,020mm
Weight	Approx. 60kg

Note: Water Supply Pressure should be 98 to 294kPa (1 to 3kg/cm²)