

CV600



Specifications

| Model | CV600 |
|-----------------------|---------------------------------------|
| Plasma mode | RIE • DP |
| Electrode | Parallel flat stage plate |
| RF generator power | 100 - 600W |
| Oscillating | 13.56MHz |
| Control system | Digital display and one touch control |
| Chamber | (W) 310 x (D) 220 x (H) 140mm |
| Electrode stage | (W) 230 x (D) 140mm |
| Chamber material | SUS |
| Reaction gas | O ₂ ,Ar |
| Purge gas | N ₂ or Dry air |
| Flow meter | Flow meter |
| Vacuum pump | Vacuum pump(120L/min) |
| Gas system | 2 line and purge line |
| External dimensions | (W) 480 x (D)500 x (H) 1,065mm |
| Power supply(50/60Hz) | AC200/220/380V Three phase 7A |