



- Rotary/Reciprocation Shaking Changeable Type for the Appropriate Mixing, Extracting and Stirring
- Various Shakings by changing Shaking Stage (Optional Item)
- Timer can adjust Shaking Time
- Shaking Frequency is Digital Display
- Shaking Stage can be changed easily
- Compact Design

MK201D with Two Layers Spring Shaking Rack

■ Specifications

Model	MK201D
Shaking Type	Rotary and Reciprocation Shaking (Manual Changing)
Shaking Pitch	Rotary / Reciprocation : 30mm for both
Shaking Frequency	20 to 200rpm, Random Adjustment, Shaking Frequency Digital Display
Timer	0.5 second - 100 hours
External Dimensions	W442 x D415 x H130mm
Weight	Approx. 19kg
Power	AC100V, 50/60Hz, Approx. 2A

■ Optional items

Shaking Rack for Erlenmeyer Flask Fixing Type



Product code	Capacity	Quantity
232170	100mL	20pcs
232171	200mL	10pcs
232172	300mL	7pcs
232173	500mL	4pcs
232174	1,000mL	2pcs

Shaking Rack for Erlenmeyer Flask Angle Adjustable, Nail Clamp, Inclination Type



For Erlenmeyer flask (1set)		
Product code	Capacity	Quantity
232175	100mL	8pcs
232176	200mL	4pcs
232177	500mL	2pcs

Two Racks Shaking is available simultaneously
(Inclination Setting Platform is necessary)

**Shaking Rack for Conical Tube
Angle Adjustable Type**



For Conical Tube (1set)
Product code 232178 Capacity 15ml for 12pcs
O.D 16.5mm x less than L110mm

Two Racks Shaking is available simultaneously
(Inclination Setting Platform is necessary)

**Shaking Rack for Centrifugal Tube
Angle Adjustable Type**



For centrifugal tube
Product code 232179 Capacity 50ml for 8pcs
O.D 30mm x less than L110mm

Two Racks Shaking is available simultaneously
(Inclination Setting Platform is necessary)

**Shaking Rack for SUS Spring
Almighty Type**



A Type: for Polyethylene Bottle
Product code 232180 Capacity 250ml max 9pcs

**Shaking Rack for SUS Spring
Almighty Type**




B Type: for Erlenmeyer Flask and Testing Tube
Product code 232181
17 x 15=255 squares (14mm square)

Sticky Non-slip sheet



Product code 232182
W400 x D330mm Material: Silicon Rubber

Two Layers Spring Shaking Rack



C-Type: for Test tube and Erlenmeyer flask
Product code 231398
Test tube 16mm dia 45deg of inclination 105pcs
Erlenmeyer flask 100ml x 16pcs 200ml x 10pcs
500ml x 4pcs 1,000ml x 2pcs

Inclination Type Setting Platform



Product code 232078
Necessary for Angle Adjustable Type
Two Inclination Type Setting Platforms are available simultaneously

*Glassware and Polyethylene Bottle in the photo are not included in the Optional Items.