

Magnetic Stirrer

M-21

Model	M-21
Revolution	200 to 2,500 rpm
Stirring Stage Size	92 x 120mm

Features

- Suitable for Titration Stirring and for Small Amount Stirring
- Super Light and Compact Type
- Stirring Stage is Aluminium with Chemical-proof



Specifications

Model	M-21
■ Standard	
Revolution	Approx. 200 to 2,500rpm
Motor Type	Shaded Pole Motor
Stirring Stage Size	W92 x D120mm
External Dimensions	W96 x D129 x H73mm
Power	AC100V, 50/60Hz, 0.5A
Weight	Approx. 0.6kg
Accessory	Stirrer Bar, 2pcs.

■ Optional items

Description	Model No.	Product Code
Stirrer Fixing Support	M-21	231198