

Magnetic Stirrer

MA100/300

Model	
100	300
Revolution	
100 to 1,500 rpm	400 to 1,550 rpm
Stirring Capacity	
50 to 1,000ml	100 to 3,000ml

Features

- **MA100**
Small Size
Loaded with DC Motor Electronic Control
- **MA300**
Long Life
Max. 3ℓ Stirring



MA100/300

Specifications

Model	MA100	MA300
■ Standard		
Stirring Stage Size*	W77 x D135mm *Not include Liquid Drop Protection Part	W178 x D165mm *Not include Liquid Drop Protection Part
Stirring Capacity	50 to 1,000ml	100 to 3,000ml
Revolution	Approx. 100 to 1,500rpm	Approx. 400 to 1,550rpm
Motor	DC Motor, Electronic Control, 0.3W	AC Motor, Shaded Pole Type, 0.7W
Power	AC100V, 50/60Hz, 0.2A	AC100V, 50/60Hz, 0.5A
External Dimensions	W84 x D150 x H60mm	W181 x D197 x H83mm
Weight	Approx. 1.1kg	Approx. 1.6kg
Accessory	Stirring Bar 30mm, 1pc.	