

Mixer for Small Container
MA300A/300B



MA300A

- Basic Type Digital Indication Magnetic Stirrer.

MA300B

- Thin Type Digital Indication Magnetic Stirrer.

Model	
300A	300B
Revolution (rpm)	
100 to 1,200 rpm	
Stirring Capacity	
50 to 3,000ml	50 to 1,000ml

■ Specifications

Model	MA300A	MA300B
Stage Size	W172 x D156mm	
Stirring Capacity	50 to 3,000ml	50 to 1,000ml
Revolution	100 to 1,200rpm	
Motor	DC Motor	
Power	AC100V, 50/60Hz, 0.5A	
External Dimensions	W172 x D156 x H53.7mm	W172 x D156 x H37.7mm
Accessory	Stirrer Bar, 30mm, 1pc.	