

Mixer for Small Container

MA300A/300B

| Model | |
|-------------------|---------------|
| 300A | 300B |
| Revolution | |
| 100 to 1,200 rpm | |
| Stirring Capacity | |
| 50 to 3,000ml | 50 to 1,000ml |

Features

- **MA300A:** Basic Type Digital Indication Magnetic Stirrer.
- **MA300B:** Thin Type Digital Indication Magnetic Stirrer.



MA300A/300B

Specifications

| Model | MA300A | MA300B |
|---------------------|-------------------------|---------------------|
| Stage Size | W165 x D157mm | W165 x D153mm |
| Stirring Capacity | 50 to 3,000ml | 50 to 1,000ml |
| Revolution | 100 to 1,200rpm | |
| Motor | DC Motor | |
| Power | AC100V, 50/60Hz, 0.5A | |
| External Dimensions | W165 x D157 x H53mm | W165 x D153 x H37mm |
| Accessory | Stirrer Bar, 30mm, 1pc. | |