

Magnetic Stirrer

MB800

Model	MB800
Revolution (rpm)	70 to 1,200
Stirring Capacity (ml)	100 to 10,000

FEATURES

- Top Plate is Chemical-proof, Anodized Aluminium Finish
- Automatic Power Shut Down Function by detecting the abnormal temperature of Water Bath



Specifications

Model	MB800
■ Standard	
Stirring Plate Material	Aluminium
Stirring Plate Dimensions	W250 x D220mm
Stirring Capacity	100 to 10,000ml
Revolution	70 to 1,200rpm
Motor	AC Motor, Electronic Control
Overheat Preventor	70 to 200 deg.C
Sensor	Themista
Service Receptacle (for Oil Bath)	AC100V, 50/60Hz, 10A
Safety Device	Electric Leakage Breaker, Overheat Preventor for Oil Bath Power Shut Down
Power	AC100V, 50/60Hz, 11A (at Operation with Oil Bath, Max.)
External Dimensions	W250 x D270 x H150mm
Weight	Approx. 4.2kg
Accessory	Stirrer Bar 40mm, 1pc.