

Magnetic Stirrer
MC801/MF820



MC801/MF820

- **MC801:**
 - With Three Rotative Direction (Counterclockwise, Reverse, Clockwise)
 - Strong Stirring with little Temps. Increasing
- **MF820:**
 - Efficient Stirring by the Stronger Magnetic Power

Model	
MC801	MF820
Revolution	
App. 80 to 1,800 rpm	App. 80 to 1,500 rpm
Stirring Capacity	
100 to 10,000ml	100 to 20,000ml

■ Specifications

Model	MC801	MF820
Performance		
Stirring Stage Material	Stainless Steel	
Stirring Stage Size	W272 x D270mm	W270 x D260mm
Stirring Capacity	100 to 10,000ml	100 to 20,000ml
Revolution	Approx. 80 to 1,800 rpm	Approx. 80 to 1,500 rpm
Motor	DC Motor Electronic Control	
Rotative direction of motor	(CCW-Counterclockwise direction), Reverse Stirring, (CW-Clockwise direction),	
Timer	Adjustable Reverse Stirring Time in between 10 to 120min.	N/A
Power	AC100V, 50/60Hz, 1A	
External Dimensions	W278 x D286 x H81mm	W276 x D295 x H80mm
Weight	Approx. 3.6kg	Approx. 3.8kg
Accessory	Stirring Bar 30mm and 40mm, 1pc. each	Stirring Bar 40mm 1pc.

■ Stirable even on a Thick Catalogue



■ Optional Item



*Glassware in the picture are not included in the Optional Item.

Description	Product Code
Spring Shaking Rack	231394