

## Magnetic Stirrers

**MG120/600**

| Model                  |                       |
|------------------------|-----------------------|
| 120                    | 600                   |
| Revolution (rpm)       |                       |
| 100 to 1,500           | 200 to 1,500          |
| Stirring Capacity (ml) |                       |
| 5 to 1,000<br>12 pcs.  | 50 to 2,000<br>6 pcs. |

### Features

- Simultaneous Stirring up to 12 samples (MG120)
- Individual Revolution Control for 6 points (MG600)
- Chemical-resistant plate with ceramic coating (MG600)



MG120/600

## Specifications

| Model             | MG120  | MG600   |
|-------------------|--|---|
| Stirrer           | 12 pcs. Interlock Type (with OFF Volume)       | With Individual Revolution Control for 6 points (With OFF Volume) |
| Plate Material    | Ceramic Coating                                |   |
| Plate Size        | W295 x D198 mm                                 | W449 x D300 mm  |
| Stirring Capacity | 5 to 1,000ml, 12 pcs. (Max.)                   | 50 to 2,000ml x 6pcs.   |
| Revolution        | Approx. 100 to 1,500rpm.                       | Approx. 200 to 1,500 rpm.   |
| Motor             | DC Motor x 12pcs. Interlock Electronic Control | DC Motor 0.7W x 6pcs.   |
| Power             | AC100V 50/60Hz Approx. 0.15A                   | AC100V 50/60Hz 0.5A   |
| External Size     | W303 x D234 x H55 mm                           | W456 x D335 x H90 mm  |
| Weight            | Approx. 4kg                                    | Approx. 7.5kg   |
| Accessory         | Stirrer Bar 25mm 12 pcs.                       | Stirrer Bar 30mm 6 pcs.   |