

## Touch Mixer

**MT-31/51**

Model	
31	51
Revolution (50/60Hz)	
2,800 to 3,300 rpm	600 to 3,000 rpm
Stirring Stage Size	
70 mm	

### Features

- Compact Type Touch Mixer for Continuous Operation
- Loaded with Stirring Stage Switch which has the mixing start by putting Sample on it.
- MT-51 is loaded with Revolution Adjustment Dial.



MT-31/MT-51

## Specifications

Model	MT-31	MT-51
<b>■ Standard</b>		
Revolution	2,800/3,300rpm (50/60Hz)	600~3,000rpm
Motor	Shaded Pole Motor, 2W	Direct Current Motor, 10W
External Dimensions	W104 x D155 x H123mm	W128 x D165 x H125mm
Weight	1.3kg	1.8kg
Switch	ON1: Touch Switch, ON2: Continuous Switch	
Stirring Ware	Testing Tube, Centrifugal Tube, Colour Comparison Tube, Erlenmeyer Flask (Max. 100ml)	
Stirring Surface Size	70mm	
Stirring Surface Material	Polyurethane	
Power	AC100V, 50/60Hz, 1A	