

## Water Purifier

**WL100**

Model	WL100
Collection Type	Ion Exchange Type
Collection Pure Water	Ion Exchange Water


### Features

- Simply Collectable Large Amount of Ion Exchange Water only by connecting to Tap Water
- Built-in Voltage Circuit in Water Quality Meter can check Water Quality immediately.
- Polyethylene made Cartridge for lesser contamination
- No AC Power is required, and Economical Price
- Easy Maintenance



This sets is WL100+Filter Stand (1)+Filter Housing Unit (2)+Active Carbon Filter+Membrane

### System Example



Possible to eliminate the other particles of Inorganic Ion

## Specifications

Model	WL100
<b>■ Performance</b>	
Collection Type	Ion Exchange Type
Collection Pure Water	Ion Exchange Water
Production Amount	Below 2 to 5L/min.
<b>■ Components</b>	
Ion Exchanger	Cartridge Type (Mixed Bed Type, Resin Amount 10L)
Water Quality Meter	0 to 10 x 10-4S/m (0 to 10μs/cm) (Analogue Display)
<b>■ Standard</b>	
Power (for Water Quality Meter)	Dry Cell Battery, DC9V
External Dimensions	Φ180 x H820mm (Not Including Coupler and Stand)
Weight (at Operating)	Approx. 18kg

### ■ Optional items


Remark	Description	Model No.	Product Code
	Connection unit A	WL100/21P	253677
	Connection unit B	WL100+Filter Housing Unit	253678
	Water supply port unit OWH10	WL21/21P, WE21, RO21/81	253686
	High purity ion exchange cartridge ICR-1	WE21,WL21/21P* *required 2 pcs.	000813
○	Stand for AS32A	WL21P	253673
○	Stand for AS32B	WL21P	253674
○	Stand for AS31A	WL21	253671
○	Stand for AS31A	WL21	253672
	Filter stand	WL100	253676
	Filter housing unit	WL100	253675
	Pretreatment cartridge PWF-1	WL300	253099
	Ion exchange resin cartridge CPC-L	WL300	253253
	Activated Carbon filter YAC-250	WL21/21P	9020026002
	Ion exchange resin cartridge ICR-2	WL21/21P/WE21	000811
	Ion exchange resin cartridge ICR-2, regeneration charge	WL21/21P/WE21	000812
	Ion exchanger, IE bomb	WL100	000821
	Ion exchanger, IE bomb, regeneration charge	WL100	000822
	Membrane filter PSECG20SSIAW1	WL21/21P/WE21	9020010007

#### Note:

The transportation cost will not be required if the order is together with main unit except the item marked "○". In case you order the item marked "○" together with the main unit, please ask Yamato to obtain the correct transportation cost.



Connection Unit G (WL100+WG262)



Filter Stand + Filter Housing Unit