

Water Purifier  
**WL100**



Model
WL100
Collection Type
Ion Exchange Type
Collection Pure Water
Ion Exchange Water


This sets is WL100+Filter Stand (1)+Filter Housing Unit (2)  
+Active Carbon Filter+Membrane

- Simply Collectable Large Amount of Ion Exchange Water only by connecting to Tap Water
- Built-in Voltage Circuit in Water Quality Meter can check Water Quality immediately.
- Polyethylene made Cartridge for lesser contamination
- No AC Power is required, and Economical Price
- Easy Maintenance

■ Specifications

Model	WL100
<b>Performance</b>	
Collection Type	Ion Exchange Type
Collection Pure Water	Ion Exchange Water
Production Amount	Below 2 to 5L/min.
<b>Components</b>	
Ion Exchanger	Cartridge Type (Mixed Bed Type, Resin Amount 10L)
Water Quality Meter	0 to 10 x 10 <sup>-4</sup> S/m (0 to 10 <sup>5</sup> µs/cm) (Analogue Display)
<b>Standard</b>	
Power (for Water Quality Meter)	Dry Cell Battery, DC9V
External Dimensions	Φ180 x H820mm (Not Including Coupler and Stand)
Weight (at Operating)	Approx. 18kg

■ System Example




Possible to eliminate the other particles of Inorganic Ion

■ Optional items

Description	Model No.	Product Code
Connection unit B	WL100+Filter Housing Unit	253678
Connection unit G	WL100+Filter Housing Unit	253668
Filter stand	WL100	253676
Filter housing unit	WL100	253675
Activated Carbon filter	WL100	9020026002
Ion exchanger, IE bomb	WL100	000821
Ion exchange resin only (without cylinder)	WL100	9110010001
Membrane filter PSECG20SSIAW1	WL100	9020010007



Connection Unit G (WL100+WG250)



Filter Stand + Filter Housing Unit